# Revolutionary device for the measurement and control of all major parameters of laboratory weathering

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Sandra Schneider - ATLAS

Atlas Material Testing Technology, a world leader in material durability testing solutions, has introduced a revolutionary device for the measurement and control of all major parameters of laboratory weathering, the All-In-One-Sensor (AIOS). This sensor mounts onto a standard Ci Series Weather-Ometer® rack, continuously and simultaneously measuring irradiance, black panel temperature, chamber air temperature, and relative humidity. Different configurations are available to measure irradiance either at 340nm, 420nm, or 300-400nm, as well as options for the reference sample temperature (uninsulated black panel (BPT) or insulated black standard (BST)). The sensor is compact, only taking the space of one standard sample holder. It is based on the popular Xenocal and Xenosensiv technologies used in the Atlas Xenotest® and Suntest® product lines.

A desirable characteristic of any laboratory weathering test is precise control of all parameters to provide better repeatability and reproducibility. Irradiance and a reference black panel temperature have always been measured at the same sample surface. However, it has been widely known that the common methodologies for measuring chamber air temperature and relative humidity vary depending on the product line and instrument manufacturer. These differences contribute to variability of test results between different types

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of instruments. Communicating the measurement and control of these parameters at the sample level reduces this variability. The irradiance sensor has a special diffuser to appropriately measure a wide field of view, including a correction for the angle of incidence to the sensor (commonly called a cosine phi correction). This correction is a requirement of ISO 9370, which is referenced in the fundamental

methods for laboratory weathering such as the ISO 4892 and ISO 16474 series of standards. To protect the optical system against the influence of the environment, a quartz dome is covering the system.

Since the sensor's electronics and filter system are temperature sensitive, all sensors are temperature compensated. An internal temperature sensor continuously measures the sensor temperature and this data is used to correct the measured data before it is transferred into a computer or the control system.

Since the sensor is exposed to the full spectrum of the xenon lamp, a temperature shield protects the sensor from excessive heating. All sensors in the system are calibrated according to national or international standards.

Irradiance is calibrated by a process traceable to a National Standards body through Atlas' optical laboratory in Germany, which is accredited by DAkkS, a nationally recognized accreditation organization in Germany. The BST and BPT sensors are calibrated by a contact method or by a new and patented contactless method using xenon lamps as source for IR radiation.

The relative humidity and CHT sensors are calibrated by the manufacturer. As a part of normal preventive maintenance, these sensors can be easily replaced on an annual basis, but Atlas is working to develop a process to calibrate these sensors as well.

The sensor is equipped with its own firmware, which controls the function of all the different sensors and external communication, as well as all the sensor calibration data, sensor identification information, and other sensor characteristics. The sensor has an internal datalogger to take offline measurements.

Therefore, the sensor could be mounted in any laboratory instrument (or outdoor). This data can then be analyzed, and test parameters could be set or modified to provide the same conditions.

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