Rapid optimisation of pigment dispersion using parallel "high throughput" grinding

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Stephen Bysouth, Ali Muhammad - AUTOMAXION

INTRODUCTION

Choice of equipment to break down aggregated pigment particles depends on factors such as pigment type, pigment particle size and distribution, pigment hardness, premix viscosity, ink vehicle, solids content in the formulation, temperature limitations, aeration tendencies and the required extent of pigment dispersion. Once primary particles are obtained, re-aggregation must be prevented usually by using additives to control surface charge (electrostatic stabilisation) or to cover the surface and 'mechanically' prevent aggregation (steric stabilisation). How the solids can best be dispersed in a medium depends on intrinsic characteristics of the pigment, etc.

For carbon black, surface oxidation improves its dispersibility but for conductive inks, any oxygen-containing surface groups, markedly reduces electrical conductivity. These multitude of factors make it difficult to disperse the pigment in the ink vehicle and for a new finished ink or coating, involves tests on a very large number of samples and a lot of time, raw materials and waste.

Over the last 15 years or so, high-throughput techniques employing hardware, software, data analysis and presentation, have been applied to formulation of materials. These have included large integrated systems but also standalone tools for specific parts of the workflow. In this article, the optimisation of conductive carbon black pigment formulations and the use of some simple tests to select suitable candidate materials and a multisample

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mill to prepare the dispersions is compared to a recirculating bead mill.

MATERIALS AND METHOD

The first phase in any HT experimental program if possible, is to reduce the number of experiments. We had the possibility of using two different, waterborne, binder-free dispersions of carbon black pigments (Tab. 1) with low volatile matter content that allows the analysis of the extent of pigment dispersion brought about principally by the action of the mills. High molecular weight, graft copolymers are recommended as dispersing additives for use with carbon black pigments possessing neutral surface. We chose to look at three examples (Tab. 2).

A simple 'sinking' test was performed to select the best candidate dispersant, where equal amounts of pigment were added to 5wt% active matter dispersant solutions in a test tube and were then left at ambient temperature without agitation. The samples were observed after 24 hrs and the dispersant solution Tego 760W had the minimum sedimentation of pigment and maximum opacity of solution

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Tab. 3 Process times for milling and grinding steps for the two mills

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and was selected for optimisation of dispersion of the pigments. The maximum amount of pigment in a dispersion that could be re-circulated properly in the re-circulating bead mill was set as the upper limit for pigment loading. However, it was decided not to proceed with Enasco 350G pigment because at high suitable loadings, it was difficult to process in either mill.

To optimise the dispersant dosage, a ladder series of dispersions was prepared using the Tego 760W for use with the Vulcan XC605 pigment.

Prior to milling the dispersions in the Mini 50 Motormill (Eiger Torrance Limited), the ingredients were premixed for 30 min at 1500 rpm using a high speed overhead stirrer. The milling chamber of the Eiger Torrance bead mill has an empty volume of 50 ml which was filled with 30-35 ml of

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1-1.25 mm zirconia milling beads.

The recommended milling time of 3 hrs and milling speed was 3750 rpm was used and we found that a minimum dispersion volume of approximately 75 ml was required to ensure proper recirculation of the dispersion in the mill. For the Automaxion Mill (Pulverisette 6 with 8 'EPA' vial adapter), no premixing was needed, saving time and material. The ingredients were loaded directly into the 20 ml vials to make a sample volume of about 10 ml. An equal volume of 1-1.25 mm zirconia milling beads were added. However, the vials were not full since the formulation occupied the empty space between the beads. The vials were processed for 4 hrs 30 mins at 550 rpm.

RESULTS

In these experiments we compared how well the multisample mill compared with the recirculating mill.

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Tab. 4 Particle size analysis results for the dispersions prepared using Automaxion mill and Eiger Torrance mill (DOWP: Dispersant On the Weight of Pigment)

Level (DOWP) of dispersant 760W with 22 wt% of Vulcan XC605 pigment <i>Livello (DOWP) del disperdente 760W con 22% di peso del pigmento Vulcan XC605</i>	Mean particle size, nm <i>Dimensione della particella principale</i>	
	Eiger Torrance mill <i>Mulino Eiger Torrance</i>	Automaxion mill <i>Mulino Automaxion</i>
12%	113.8	90.47
15%	119.3	90.26
18%	123.3	90.23

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The process step times (Tab. 3) were typical for these mills and the milling/grinding times for each mill were developed to provide the required properties of the pigment concentrates created.

If single samples are being made, then there is no advantage in using the Automaxion mill since it requires more time for the actual milling and exceeds the extra time required to clean the Eiger Torrance mill.

It could be argued that the Automaxion mill does generate vials that need to be cleaned, but we assumed zero time for that, since that cleaning can be done while the mill is running other samples, unlike the Eiger Torrance mill.

However, it is obvious that the Eiger Torrance mill is limited to one sample per day (two at a push), whereas the Automaxion mill can comfortably process 8 samples a day.

The volume of sample made should be the minimum possible but sufficient for the needs of testing. In our case, we needed samples for particle sizing by the Malvern Zetasizer (a few drops), rheology using a TA Instruments AR-1500 EX rheometer (a few mls) and stability testing via storage and repeated rheology (up to 10 mls).

The Eiger Torrance mill produced about 75 ml per sample,

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Tab. 5 Zeta potential results for the dispersions prepared using Automaxion mill and Eiger Torrance mill (DOWP: Dispersant On the Weight of Pigment)

Level (DOWP) of dispersant 760W with 22 wt% of Vulcan XC605 pigment <i>Livello (DOWP) del disperdente 760W con 22% di peso del pigmento Vulcan XC605</i>	Zeta potential (mV) <i>Potenza zeta (mV)</i>	
	Eiger Torrance mill <i>Mulino Eiger Torrance</i>	Automaxion mill <i>Mulino Automaxion</i>
12%	-23.2	-17
15%	-22.3	-29.5
18%	-21	-33.4

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which is of course more than enough for these tests.

The Automaxion mill only produced about 10 ml of sample which was also enough but it is also possible to rapidly produce duplicate samples if needed.

The throughput of the mill is irrelevant if the milling/grinding results are not comparable. In most practical situations there are differences in performance between laboratory mills or against a production mill but these are known and well understood.

For these two mills, each have variable parameters such as milling speed, time, media size, media fill or pumping rate for the Eiger Torrance mill. As yet we have not invested time in finding how all the parameters for the two mills correlate but the stated conditions were the result of some initial experiments and experience.

For these inks, the important metrics are the particle size distribution, zeta potential, the rheology and the rheological stability after a period of storage. For particle size distribution we used the Zetasizer Nano ZS (Malvern Instruments). 0.25

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g of each sample was diluted with 50 ml of de-ionised water and sonicated for 1 min. The refractive index and absorption of the pigment were set at 1.675 and 1.00, respectively.

Due to the extent of dilution, the default values of 1.33 and 78.5 were used for refractive index and dielectric constant of the dispersant (i.e., water), respectively.

The sample was equilibrated at 25°C for 60 sec in the disposable sizing cuvette (DTS0012) prior to analysis.

The zeta potential of an aliquot of the same diluted sample was determined by performing an electrophoresis experiment in a disposable zeta cell (DTS1060C).

The particle size data in Tab. 4 shows that the Automaxion mill brought about a greater decrease in the mean particle size and was consistent for each sample, compared to the Eiger mill.

This is confirmed by the distribution plots (Fig. 1). This decrease may be due to the slightly longer milling duration in the case of the Automaxion mill or, more likely, the greater effectiveness of the planetary milling action of the Automaxion mill which is also confirmed as there appeared to be some residual poorly ground material in some of the Eiger Torrance mill samples (Fig. 1), which skews these results. The results of zeta potential measurements (Tab. 5), show

Fig. 1 Particle size distribution plots for 22% VXC605 pigment loading and the three level of dispersant 760W measured soon after processing

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Fig. 2 Viscosity profiles of dispersions of 22% VXC605 pigment with (a) 12%DOWP 760W, (b) 15%DOWP 760W and 18%DOWP 760W prepared on Automaxion and Eiger Torrance mills measured soon after processing

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Fig. 3 Viscosity profiles of dispersions of 22% VXC605 pigment with (a) 12%DOWP 760W, (b) 15%DOWP 760W and 18%DOWP 760W prepared on Automaxion and Eiger Torrance mills, measured soon after processing and after 4 weeks

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considerable difference between the series of dispersions processed on the two mills. This may be a result of the consistent particle size and, since the particle size is consistent for the Automaxion mill and the resulting zeta potentials may better show the effect of varying dispersant levels.

However, since the Dispersant 760W is known to sterically stabilise pigment particles, further studies would be needed to elucidate the relationship further. Of course,

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these samples had to be diluted for measurement which is known to alter the characteristics of the test sample, whereas samples for rheological analyses were carried out on actual dispersions. For rheometry, we used a TA Instruments AR-1500 EX cone and plate rheometer with 60 mm stainless steel cone with a 1° taper angle and 27 microns truncation gap. Prior to measurement, each sample was manually shaken for 2 min and then left at ambient temperature (22-25°C) for 5 min. Measurements were made without pre-shearing the sample. A continuous ramp test was employed and the viscosity profile was recorded over a shear rate range of 1/s to 1000/s.

The results of the rheology analysis Fig. 2 correspond with the particle size and zeta-potential results in that the viscosity profiles for the Automaxion mill are consistent whilst those for the Eiger Torrance mill show greater variation, as do the particle sizing results. For inks, it is also important to understand the stability of the dispersions, since post-milling build-up of viscosity is possible if the pigment milling is not done effectively and optimally or the dispersant is too much in deficit or excess. For stability analysis, the dispersions were stored at ambient temperature (22-25 °C) and the viscosity profiles were recorded after four weeks. These viscosity profiles are presented in Fig. 3. The viscosity of each of the three dispersions was found to be very stable. This showed that the Automaxion mill produced very stable dispersions of Vulcan XC605. If the viscosity profiles were examined on their own, it might be concluded that the level of dispersant was equally effective between 12 and 18 % when looking at the results using the Eiger Torrance mill alone. However, if the consistent particle size produced by the Automaxion mill are considered with the viscosity stability results, there is some evidence that 15% of the dispersant produces more stable formulations. This is masked by the variation in particle sizes obtained using the Eiger Torrance mill.

CONCLUSIONS

Milling by the Automaxion mill is comparable to the re-circulating bead mill and may be better since, in a recirculating mill a portion of the formulation being processed is always in the recirculating loop which can lead to broader distributions and unmilled sample. The Automaxion system offers several advantages over a recirculating mill including:

1. Less raw materials required and consequently, less waste.
2. Less time because several samples can be processed at the same time.
3. Less cleaning time and milling directly in the vials in which the dispersion can be stored.

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