Dispersion system for laboratory experimentation to grind or disperse powders in suspension

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Dr. S. Bysouth
Automaxion

Automaxion Sarl based in France, has developed a bead milling and dispersion system for laboratory experimentation that can grind or disperse powders in suspension. Originally developed for comminution of crop protectants in a sophisticated robotic system, the company has been broadening the applications with its customers to include pharmaceuticals, catalysts and, of interest to the European Coatings readers, pigments for coatings and inks.

Based on a patented variation of planetary milling, the system was developed to grind multiple samples in vials which could be handled by a robot but when used unautomated, as most users do, the vials and the way they are held in the system still present considerable advantages when processing pigments. Sample size can be varied by selecting adapters for vials between 20 and 100 ml (Fig. 1). Below 20 ml, it is more harder to get energy into the beads (though for pharmaceutical applications, 2 ml vials have proven useful) and above 100 ml, a more conventional planetary mill is preferred as the vials cannot be accommodated in the special adapters.

When considering performance, the particle size distribution is comparable or superior to such things as recirculating or basket mills. Of course, this is dependent

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Fig. 1
Sample size can be varied by selecting adapters for vials between 20 and 100 ml

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Fig. 2
Bead milling and dispersion system for laboratory experimentation that can grind or disperse powders in suspension

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on the way the system is operated but with the right beads and conditions, narrow distributions can be obtained with a D_{50} of 0.1 μm or less (which is the practical limit for bead milling devices (fig. 2)). In some applications, this can be achieved in a shorter time compared to the conventional technique. For example, in a solvent coatings application, a milling time of about 10 minutes has been used compared to about 3 hours using a

laboratory basket mill. This application highlights another advantage: since the milling is carried out closed vials, there is no evaporation during milling which means that at the end of the milling the formulation is consistent with the original recipe, whereas evaporation from an open basket mill over 3 hours means that the solvent content of the final dispersion is unknown. In fact, the mill can handle pressure vessels and has been used for milling in a highly volatile solvent. Of course for this application and in any volatile solvent, the vial must be capable of withstanding the pressure at the temperatures generated by the milling. However, heat build-up is also limited by the way the vials are held in a large heat-absorbing and dissipating aluminium block, the comparatively small sample size and the air circulation around the milling space and adapters. Furthermore if temperature is critical, the mill can be operated with fewer samples, at a reduced speed or programmed to allow rest periods but as yet, in our experience, this has not been needed.

For aqueous systems and with good quality vials and closures, our experience is that vial leakage or breakage is very unusual. For volatile solvent systems, again as long as the vial chosen is suitable, leakage or breakage is rare but, especially if the solvent is flammable, secondary containment can be provided using an O-ring on the vial adapter. One word of caution, to avoid breakage, the vials used must be those originally specified for the adapters.

Saving the best news 'till last, the biggest advantage with the system comes from the sample throughput, which is especially important when working on understanding formulation variables, optimising ingredient concentrations and troubleshooting. Of course, the number of samples that can be processed at the same time, depends on the model of the mill but for a 4 adapter position mill and 30 ml vials, 32 samples can be milled at the same time. In addition, since samples are held within vials, there is nothing to clean and vials can be exchanged rapidly. In addition, you don't need to do any predispersion outside the mill since this can be carried out (if desired) by initially running the mill at a reduced speed, but we have found this step to rarely be needed. In one study

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where the University of Leeds was dispersing carbon black for special inks, we found that compared to a recirculating bead mill, the Automaxion system would allow eight samples to be processed per day as opposed to only one. Because of the small sample size, it is best to use critical tests on these samples first. If you then find you need a larger volume, you can use the mill to make replicates to combine for such tests. Of course, if you produced 32 samples at conventional volumes, you would also consume considerable quantities of ingredients and equivalent waste. Recently, we worked with a customer who uses valuable pigments and small sample size was a critical factor in using the Automaxion mill.

As with any new piece of equipment, there is some learning, adaptation and method development needed when first using the system but of course, since the mill rapidly produces samples, even this is relatively fast. Of course, Automaxion can help and if you want to try the mill on your samples, we can run some for you (you send everything to us already in the vials and we run them, returning them to you unopened) or you can have a mill in your laboratory for a period for testing. In either case, we support you through a two-way discussion to trouble shoot and help you understand how to use the system to best advantage. Alternatively, we can come to your laboratory to assist.

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