Analytical measuring instruments characterize particles and dispersions

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Arnold Uhl - LUM

The product portfolio of Berlin based LUM GmbH provides innovative analytical instruments for the characterization of the entire life cycle of materials, which help formulators to manufacture paints, varnishes, adhesives and building chemicals. For the quality control of raw materials, on the basis of their particle and or droplet size distributions, the dispersion analyser LUMiSizer is used for nanoparticles and the separation analyser LUMiReader PSA is used for microparticles. For the characterization of dispersibility and separation stability or shelf life the user needs to know the product properties in its original concentration. This requirement is fulfilled by both instruments, too. Depending on the formulation properties and application the samples are measured in real-time or directly physically accelerated, to obtain the information in time lapse mode, instead of waiting for weeks, months or even years. The LUMiSizer is an evolutionary development based on the LUMiFuge platform with some significant additional capabilities that instantaneously measures the extinction (space- and time-resolved) of the transmitted light across the entire length of your sample using the STEP-Technology. Because of an enhanced optical system, the instrument lets you analyse particle and droplet velocity distributions for creaming and sedimentation phenomena without the need of any material data and performs particle sizing (ISO 13318- 2).

Run 12 samples at a time under a wide range of viscosities, temperatures, and real-world concentrations. Product is the only instrument in the world that allows you to obtain particle velocity distributions without having to know any material constants. For paints and varnishes industry analyse aqueous, non-aqueous, Newtonian or non-Newtonian systems – the choice is yours. You can measure the effects of concentration, shape, and

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colloidal forces on your particle size distribution and stability. You can create your own hindrance functions or work with published ones. Either way, you can accurately measure and predict ideal and non-ideal particle behaviour. The exceptional precision of the system ensures that all particle size fractions are measured in a short time. Immediate data collection across the entire sample ensures exceptional results, whereby you can measure wide particle size distributions, determine non-idealities and gain the ability to resolve multimodalities. Applying Centrifugal Adhesion Testing (CAT)-Technology the LUMiFrac instrument determines adhesives systems and bonding strengths of composites. So it is used not only for coatings and treated surfaces but also for the characterization of composite materials, e.g. in automotive or aviation industries. Featuring a large force range and a clamping-free sample preparation, the LUMiFrac offers a wide application range. The simultaneous analysis of up to eight samples provides the operator with a high sample throughput under identical test conditions.

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