# UV curing under control with Photo-DSC

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### INTRODUCTION

The main advantages of the curing process by UV (ultra violet) light are evident: fast curing of solvent-free resins at moderate temperatures. This article presents capable thermoanalytical instruments, which are helpful tools for investigating the UV curing process and the viscoelastic properties of cured resins or of bonded parts. The following example shows the practical application example demonstrating the high impact of Photo-DSC (Differential Scanning Calorimetry) employed in the fields of research & development, quality assurance and failure analysis.

### ADVANCED PHOTO-DSC FOR UV CURING

The extension of a Differential Scanning Calorimeter (DSC) with a high-pressure mercury lamp for UV curing are known for several years. However, the quality and reproducibility of the measurement results are often poor due to the limited capabilities of the measuring device. Additionally, the competitive reaction with oxygen influences the pure UV curing. Finally, dual-cured adhesives are influenced by a reliable intensity of UV light at an accurately controlled temperature.

NETZSCH Analyzing & Testing offers

a gas-tight DSC system, which can be combined with nearly any commercial lamp via two glass-fiber guide-lines for the sample and reference (Fig. 1). The irradiation time is triggered and controlled via the measurement software. The adjustable automatic cover lift always ensures equal distance of the UV-light outlet and the sample surface. The use of an automatic sample changer (ASC) provides a high sample throughput for these fast measurements. With the application of a DSC curve subtraction and time offset, the pure resulting reaction enthalpy for UV curing can be reliably evaluated.

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Table 1 presents the main technical features of the Photo-DSC 204 F1 Phoenix®. It should be pointed out that the Photo-DSC can be retrofitted to a standard DSC at any time if the wide temperature range from -180 to 700°C is needed.

### UV CURING OF A FREE RADICAL ADHESIVE BY PHOTO-DSC

UV curing of a free radical adhesive based on acrylate was investigated by means of the Photo-DSC. The average sample input amounted to 10.51 +/- 0.08 mg. Open aluminum pans were used in an inert gas atmosphere (nitrogen). UV pulses with a light intensity of 1 W/cm<sup>2</sup> at irradiation times of 1 s were applied on the adhesive measured at different isothermal temperatures: 35°C, 45°C and 55°C.

Figure 2 illustrates the influence of temperature on the UV curing of the acrylate adhesive. In Table 2, the measured residual enthalpies with the progress of curing are listed. After an irradiation time of totally 10s, the curing enthalpy shows the highest value of 238 J/g at an isothermal temperature of 55°C. This should yield the highest network density and therefore, the highest tightness, hardness and stiffness. The total enthalpies of the acrylate adhesive cured at 45°C (230 J/g) and 35°C (228 J/g) are lower. These moderate temperatures should provide a lower tightness but a higher toughness and flexibility, which are often the target for bonding of technical parts. In practice, the same acrylate type can be tailored for different applications by applying different curing temperatures at equal UV parameters.


Furthermore, the irradiation time was modified in order to investigate the reactivity of the UV curing. The acrylate adhesive was irradiated for 2 min at the same different temperatures of 35°C, 45°C and 55°C. The average sample mass here amounted to 10.38 +/- 0.03 mg. The other parameters were kept constant.

The start of UV curing of the acrylate adhesive is shown in Figure 3 with the DSC curves (solid lines) at the three different isothermal temperatures. Here, the first derivatives of each DSC curve (DDSC, dotted lines) are also displayed.

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Their peak times and peak values represent the reactivity of the UV curing reaction. The fastest time of 6.2 s with the minimum peak value of -508 mW/mg/min provides the highest reactivity for the acrylate adhesive at the curing temperature of 55°C (red curves), followed by 45°C (6.3 s, -455 mW/mg/min, green curves) and by 35°C (6.4 s, -388 mW/mg/min, blue curves).

The conclusion of this result is that UV curing of this free radical acrylate adhesive is definitely depending on heat and, not only on UV irradiation time and intensity. Therefore, we can call this a dual cure reaction. Finally, the way of the UV irradiation process itself was investigated at a curing temperature of 35°C. The sample masses amounted to 10.427 mg (for an irradiation time of 5 x 1 s) and 10.428 mg (for 1 x 5 s).

The other parameters were the same as described earlier. The main question can easily be answered: there is no significant difference whether the acrylate adhesive was cured with one pulse of 5 s or with five pulses of 1 s each. Table 3 summarizes the equal enthalpy values for the different irradiation processes after total time of 5 s at 35°C.

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
#### CONCLUDING COMMENTS

Photo-DSC is a useful tool for the investigation of UV curing of free radical adhesives. Besides the determination of the influence of temperature on the UV curing, it also allows for the observation of the network density and therefore, the tightness, hardness and stiffness.

The influence of temperature on the UV curing process can be clearly observed. From this observation, one could draw conclusions on the hardness and stiffness of the material, especially in combination with investigations by means of DMA (dynamic mechanical analysis).

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